




# Intro to RNA-Seq

Dr. Princess Rodriguez


**MMG232**  
**Spring 2023**



normal


You have been given cells to culture:


Normal cells need their media changed every 3 days




mutation in  
*Kras*<sup>G12C</sup>

Other cells need their media changed every day!







normal



Normal + **G12C inhibitor**



mutation in  
*Kras*<sup>G12C</sup>



mutation in *Kras*<sup>G12C</sup>  
+ **G12C inhibitor**

# Typical RNA-Seq Vignette

You and your PI decide to uncover the potential mechanism of the **G12C inhibitor** and understand if there are any off-target effects.


The goal is to take the G12C inhibitor to animal studies to see if it works the similarly *in vivo*!




# What is RNA-Seq?

- ❖ Technique used to explore and/or quantify gene expression within or between conditions of an organism


## Normal Cells




## Mutated Cells






Each cell has a bunch of  
chromosomes



**Normal Cells**




**Mutated Cells**





Each chromosome has  
a bunch of genes




**Normal Cells**




**Mutated Cells**






Some genes are active more than others




## Normal Cells




## Mutated Cells





## Normal Cells





## Mutated Cells



RNA-Seq tells us which genes are active, and how much they are transcribed!





## Normal Cells




So, we can use RNA-Seq to  
measure gene expression in  
normal cells ...

## Mutated Cells






... then use it to measure  
gene expression in  
mutated cells

## Normal Cells




Then we can  
compare the two  
cell types to figure  
out what is  
different in the  
mutated cells!



## Mutated Cells



## Normal Cells



## Mutated Cells




Differences apparent for Gene 2 and  
to a lesser extent Gene 3

# Other common uses of RNA-Seq

- Expression variation in response to environmental stimuli
- Which genes are expressed in which tissues & how much?
- Study of genotype  $\leftrightarrow$  phenotype mapping
- Discovery/annotation of genes & transcripts
- Gene regulatory networks (co-regulated genes)
- Medical diagnostic

# Transcriptomics has led to the identification of new subtypes in ALL



Gu et al., 2019



# **There are 3 Basic Steps in performing RNA-Seq yourself**

1) Prepare a sequencing library

2) Sequence


3) Data analysis

# RNA-Seq experimental workflow


1

Samples of interest

A




PMID: 27548618




2

Isolate RNAs




3



Library build



5




Reads (R1 and R2)  
generated




Illumina sequencing  
versus  
ONP

4



# RNA-Seq processing



# What is a FASTQ?

```
@NS500177:196:HFTTTAFXX:1:11101:10916:1458 2:N:0:CGCGGCTG
ACACGACGATGAGGTGACAGTCACGGAGGATAAGATCAATGCCCTCATTAAGCAGCCGGTGTAA
+
AAAAAEEEEEEEEEEEEEE//AEEEEEEEEEEEEEEEEEE/EE/<<EE/AEEEAEE///EEEEAEAEAE<
```

1

2

3

4

**Each sequencing “read” consists of 4 lines of data :**

- ① The first line (which always starts with ‘@’) is a unique ID for the sequence that follows
- ② The second line contains the bases called for the sequenced fragment
- ③ The third line is always a “+” character
- ④ The fourth line contains the quality scores for each base in the sequenced fragment